

# LARYVOX® DUO BRUSH & LARYVOX® STOMA LIGHT

- ▶ For gentle and effective daily cleaning of voice prostheses
- ▶ With a brush to clean the channel of the voice prosthesis in the tracheoesophageal fistula and a sponge swab to remove deposits and residual secretions from the tracheal flange of the voice prosthesis
- ▶ With LARYVOX® Clip to allow the fixation of a stoma lamp (LARYVOX® Stoma Light, available separately). The combination of brush and lamp in one tool allows effective cleaning with one hand only and frees the other hand for other operations
- ▶ With integrated mirror on the backside of the box to release the patient from carrying extra equipment for daily care
- ▶ LARYVOX® Stoma Light, battery driven (battery include in scope of delivery), with bright LED-light to facilitate effective and safe cleaning of the voice prosthesis



## Order information:

**LARYVOX® Duo Brush**, Box with 6 pieces, Size 6 - 18

When ordering, please state the size after the item number!

**LARYVOX® Stoma light**

**REF 29110**

**REF 69101**

- Keep away from sunlight
- Keep dry
- Single patient use
- Consult instructions for use
- CE marking

Size LARYVOX® Duo Brush	Size voice prostheses	REF
6	4 - 6	29110-06
8	7 - 8	29110-08
10	9 - 10	29110-10
12	11 - 13	29110-12
18	14 - 18	29110-18